DE	<b>Chirurgische Instrumente</b> Wiederaufbereitungsanweisung	Seite 2
EN	<b>Surgical Instruments</b> Reprocessing Instructions	Page 14
FR	<b>Instruments chirurgicaux</b> Instructions de retraitement	Page 26
IT	<b>Strumenti chirurgici</b> Istruzioni per il ricondizionamento	Pagina 38
ES	<b>Instrumental quirúrgico</b> Instrucciones de reutilización	Página 50
PT	<b>Instrumentos cirúrgicos</b> Instruções de reprocessamento	Página 62
EL	<b>Χειρουργικά εργαλεία</b> Οδηγίες επανεπεξεργασίας	Σελίδα 74
SK	<b>Chirurgické nástroje</b> Návod na renovovanie	Strana 86
CZ	<b>Chirurgické nástroje</b> Pokyny na opätovné spracovanie	Strana 98
NL	<b>Chirurgische instrumenten</b> Instructies voor herbewerking	Pagina 110

# DEUTSCH

## 1. ZU DIESEM DOKUMENT

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#### ACHTUNG

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#### ⚠ VORSICHT

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### 1.1 Depictions in the document

#### ATTENTION

Warns of the possibility of material damage.

#### CAUTION

Warns of the possibility of minor personal injury.

#### WARNING

Warns of the possibility of serious or fatal injury.

#### DANGER

Warns of the imminent danger of serious or fatal injury.

#### CAUTION

This is a warning symbol which warns of dangers that could arise throughout the entire subsequent process.

– In this case, the warning symbol indicates possible consequences.

→ This is a precautionary measure against the danger.

## 2. FIELD OF APPLICATION www.biedermann.com/eifu

These instructions for use are valid for all Class I and Ir instruments from Biedermann Motech GmbH & CO.

The surgical instruments are delivered non-sterile. Before first use and before any further use, the surgical instruments must be cleaned, disinfected, and sterilized in accordance with the instructions in this document.

## 3. INTENDED PURPOSE

The Biedermann Motech instruments are divided into the following groups and intended purposes:

### Group 1: Instruments for preparation and assembly

These instruments are used for preparation and assembly of the implants, to prepare the implants and customise them, if needed, so that they can be implanted. This includes the assembly of the implants, among other things.

### Group 2: Instruments for access and preparation

The instruments are used for preparation of the surgical area and to prepare the vertebral body, in order to be able to insert the implant system. This includes, among other things, the clearing of the spinal disc space and the drilling and tapping into the vertebral bodies.

### **Group 3: Instruments for insertion, positioning, and fixing**

The instruments are used for implantation through insertion, positioning and fixing of the implants. With the help of the instrument set, the implants are inserted into the body, positioned within, and fixed at the end. The instruments include screwdrivers, head inserters, etc.

### **Group 4: Reprocessing and storage of the implant systems**

The instruments are used for storage of the implant systems and, if needed, they work as a sterile barrier in order to preserve sterility. These include, for example, caddies, sieves and sterilisation containers.

The instruments may only be used for their intended purpose in the medical field, by properly trained and qualified personnel. The treating medical doctor or user is responsible for selection of the instrument set for specific applications or surgical use, for appropriate training and information, and for adequate experience in the handling of the instrument set. Surgical instruments are intended to be standard instruments for surgical interventions in general surgery.

For this reason, we assume that the relevant laws, norms and recommendations (e.g. from the RKI or the AKI) are known.

## **4. DURATION OF APPLICATION**

Under normal conditions, the instruments are intended for an uninterrupted use for a period of less than 60 minutes. The period of application is therefore defined as temporary in accordance with paragraph 1.1, Annex IX of Directive 93/42/EEC and Paragraph 1.1 of Chapter I of Annex VIII of the European Medical Devices Regulation ((EU) 2017(745)).

## **5. COMBINATION WITH OTHER PRODUCTS**

If instruments are put together again after their disassembling, the separate parts must not be replaced by other manufacturers parts, also if a part is exchangeable because of the product's specific function (e.g. different work inputs).

## **6. WARNINGS**

The instruments are generally delivered non-sterile. After the products receipt, please check their identity, completeness, intactness and functioning before you reprocess them. It's necessary to inspect the instruments for breakage, deformation, damage and proper operation before every input. Thereby you particularly have to check parts such as cutting edges as well as all movable parts. Instruments that are fully worn, corroded, deformed, porous or otherwise destroyed have to be rejected.

## 7. REPROCESSING

Please also observe the valid legal provisions in your country and the hygiene regulations of the medical practice or hospital. This applies in particular to the various requirements as regards effective prion inactivation (not relevant for USA).

### WARNING

Procedures for prion inactivation lacking.

- With intended use, contact with patient areas relevant in terms of prions cannot be entirely ruled out in all cases.
  - High risk to the patient due to insufficient prion inactivation.
- Observe the applicable national requirements. These requirements may also include procedures which go beyond the present reprocessing instructions.

### 7.1 Instruments - Cleaning and Disinfection

#### 7.1.1 Pre-treatment directly after use

### WARNING

Residual contamination on contact surfaces.

- Patients can become infected following contact with contaminated instruments.
  - Residual cleaning agents can damage instruments.
- Disassemble instruments for pre-treatment as far as possible.  
→ Make sure that no soiling is allowed to dry onto instruments.

#### 7.1.2 Pre-treatment in the Central Sterile Services Department

### WARNING

Residual contamination on contact surfaces.

- Patients can become infected following contact with contaminated instruments.
  - Residual cleaning agents can damage instruments.
- Disassemble instruments for pre-treatment as far as possible.  
→ Make sure that no soiling is allowed to dry onto instruments.

#### Pre-treatment procedure

Selection criteria for the operating resources for the pre-treatment procedure, see chapter 7.6.1 „Operating resources for pre-treatment of instruments“.

1. Disassemble the instruments as far as possible.
2. Make up a fresh cleaning solution in the ultrasonic bath.  
When doing so, observe the requirements of the cleaning agent manufacturer, e.g. temperature, concentration, exposure time and final rinsing.
3. Take instruments out of the transport trough and manually remove coarse soiling from the instruments within max. 2 hours. When doing so, use a fresh single-use cloth soaked in cleaning solution.
4. Place the instruments in the ultrasonic bath so that they are covered with cleaning solution. The instruments should not come into contact with each other. Treat the instruments in the ultrasonic bath as follows:
  - Carefully brush the instruments.

- Rinse out the instruments. Rinse out the instruments according to the relevant instructions in the chapter 7.1.8 „Instructions for instrument groups“. For group allocation of the instruments, see the surgical instructions relating to the implant system.
  - At the same time, manipulate the moving parts and joints at least 3 times.
  - In the case of instruments with catch mechanisms, disengage where possible.
5. Activate the ultrasound for the exposure time specified by the manufacturer.
  6. Rinse off the instruments under running water at least 3 times for at least 1 minute, turning several times.
  7. Deposit the instruments in a cleaning tray.

### 7.1.3 Manual cleaning and disinfection

#### WARNING

Residual contamination due to inadequate cleaning and disinfection.

- Patients can become infected following contact with inadequately cleaned and disinfected instruments.
- Do not use manual cleaning and manual disinfection.

### 7.1.4 Machine cleaning and disinfection

#### WARNING

Residual contamination on contact surfaces.

- Patients can become infected following contact with contaminated instruments.
  - Residual cleaning agents can damage instruments.
- Disassemble instruments as far as possible.
- Trials and caddies must be machine-cleaned and disinfected separately from each other.

#### ATTENTION

Damage and contamination of instruments due to the use of an unsuitable disinfection method.

- Do not perform chemical disinfection.

### Machine cleaning and disinfection process

Selection criteria for the operating resources for the machine cleaning and disinfection process, see chapter 7.6.2 „Operating resources for machine cleaning and disinfection of instruments“.

- Place instruments in the cleaning tray with the opening facing down diagonally (for further details see table in chapter 7.1.8 „Instructions for instrument groups“).
- Place the cleaning tray in the cleaning and disinfecting device.
- Check that the instruments are not in contact with each other.
- Start the program.
- Remove the instruments after the end of the program.

### 7.1.5 Inspection

Visually inspect the instruments for the following defects:

- Corrosion
- Surface damage
- Chipping
- Soiling: Clean and disinfect soiled instruments again.

#### **WARNING**

Impermissible reuse of damaged and soiled instruments.

- Patients can become infected following contact with contaminated instruments.
  - Damaged instruments may fail.
- Send back damaged instruments or dispose of them.
- Clean and disinfect soiled instruments again or send them back or dispose of them.

### 7.1.6 Care of the instruments

#### For the European market:

Selection criteria for the operating resources for the oiling, see chapter 7.6.3 „Operating resources for the care of instruments“. Grease all threaded surfaces and joints of the instruments sufficiently with instrument oil (based on white paraffin oil) to maintain their functionality.

#### For the US market:

Optional lubrication of the instruments can ensure their long-term functionality. This is at the discretion of the user. Selection criteria for the operating resources see chapter 7.6.3 „Operating resources for the care of instruments“.

### 7.1.7 Packaging

- Assemble instruments which need to be assembled for sterilization (see associated set images and surgical instructions).
- Carry out a functional check for the instruments.
- Place the products in the sterilization tray.
- Place the sterilization tray in the sterilization container.
- When using sterilization containers, comply with national regulations; they include EN ISO / ANSI/AAMI/ISO 11607 / EN 868-8.
- For containers without integrated filter system:  
The sterilization container must be wrapped in two layer of 1-ply polypropylene wrap using sequential envelope folding techniques.
- or
- Place sterilization container into rigid sterilization container with integrated filter system and close lid (container in container system).

For requirements on sterilization containers, see chapter 7.6.4 „Operating resources for sterilization“.

### 7.1.8 Instructions for instrument groups

<b>Instrument groups</b> (for classification, see system-specific surgical instructions)	<b>Pre-treatment</b>	<b>Machine cleaning and disinfection</b>
<b>Standard</b>	see Chapter 7.1.2	Standard procedure
<b>Group A</b>	Thoroughly clean the lumen / cannulation with a brush at least 3 times	Place above a suitable rinsing mandrel
<b>Group B</b>	Brush carefully at least 3 times around the joint	Place above a suitable rinsing mandrel
<b>Group C</b>	Rinse out cannulation with syringe min. 5 ml and cannula with $\varnothing$ 1.2 mm at least 3 times	Standard procedure
<b>Group D</b>	Brush out cannulation at least 3 times and rinse out with syringe min. 5 ml at least 3 times	Place above a suitable rinsing mandrel
<b>Group E</b>	Brush out cannulation at least 3 times and rinse out with syringe min. 40 ml at least 3 times	Place above a suitable rinsing mandrel
<b>Group F</b>	Rinse out blind lumen actively by immersing at least 3 times	Standard procedure
<b>Group G</b>	Rinse out actively at least 3 times by immersing Rinse blind lumen / variable lumen with syringe min. 5 ml at least 3 times	Standard procedure
<b>Group H</b>	Brush variable lumen at least 3 times Rinse out actively at least 3 times by immersing	Standard procedure
<b>Group I</b>	Standard procedure Remove lid	Standard procedure
<b>Group S1</b>	Rinse out lumen / blind lumen with syringe min. 5 ml at least 3 times	Standard procedure
<b>Group S2</b>	Rinse out blind lumen with syringe min. 2 ml and cannula with $\varnothing$ 1.2 mm at least 3 times	Standard procedure
<b>Group S3</b>	Brush out cannulation at least 3 times and rinse out with syringe min. 10 ml at least 3 times	Place above a suitable rinsing mandrel
<b>Group S4</b>	Rinse out cannulation with syringe min. 5 ml and cannula with $\varnothing$ 0.9 mm at least 3 times	Standard procedure
<b>Group S5</b>	Rinse out annular gap with a syringe min. 2 ml at least 3 times	Standard procedure
<b>Group S6</b>	No brush Rinse actively by immersing at least 3 times. When rinsing, bend open the tight gaps at least 3 times and also rinse inner lumen	Standard procedure
<b>Group S7</b>	No brush Manipulate the moving parts backwards and forwards at least 5 times Rinse out actively at least 3 times by immersing	Positioning in fully extended state Place above a suitable rinsing mandrel

## 7.2 Sterilization

### WARNING

Use of impermissible sterilization procedures.

- Damage to implants and instruments.
- Do not use the following sterilization procedures:
  - Flash sterilization procedures
  - Gravity-displacement method
  - Hot air sterilization
  - Radiation sterilization
  - Formaldehyde sterilization
  - Ethylene oxide sterilization
  - Plasma sterilization

### WARNING

If maximum weight for sterilizing unit is exceeded, as per DIN EN 285 defined (Weight indication without sterilization container):

- Implants and instruments will not be sufficiently sterile.
- Do not exceed the following weights for container **with** integrated filter systems:
  - 1/1 sterilizing unit max. 10 kg
  - 3/4 sterilizing unit max. 7 kg
  - 1/2 sterilizing unit max. 5 kg

Note: For the container in container systems, the inner container (container without integrated filter system) is included in the sterilizing unit.

- Loading of sterilization containers **without** integrated filter system according to the given marking.
- If no marking is available, the products must be stored in one layer without mutual contact.

### Sterilization procedure

Selection criteria for the operating resources for the sterilization procedure and requirements on the sterilizer, see chapter 7.6.4 „Operating resources for sterilization“.

- Place the sterilization container in the sterilizer.
- Steam sterilization as fractionated vacuum method with at least three vacuum cycles and adequate product drying.

Recommended validated sterilization parameters:

Time	Temperature	Drying time
4 minutes	At least 132 ° (270 °F)	At least 40 minutes
3 minutes	At least 134 ° (273 °F)	At least 40 minutes

### 7.3 Storage

Store implants and instruments in sterile packaging under the following ambient conditions:

- dry
- dust-free
- protected from damage

### 7.4 Reusability

- Implants can be reprocessed, as long as they are undamaged and not contaminated, and due appropriate care is taken.
- Reuse of the instruments is envisaged. To maintain the relevant functions of the instruments, they must be treated with care and handled according to the instructions.

### 7.5 Material resilience

#### **WARNING**

Damage and contamination of implants and instruments due to unsuitable cleaning and disinfecting agents.

- Residues can result in risks for the patient.
  - Damaged implants and instruments may fail.
- Use neutral, enzymatic or slightly alkaline cleaning agents
- Cleaning and disinfecting agents which contain the following substances must NOT be used:
- Organic, mineral or oxidising acids with a minimum permissible pH value of 5.5
  - Strong alkalis, a pH value of max. 11 is permissible
  - Oxidisers (e.g. hydrogen peroxide)
  - Halogens (e.g. chlorine, iodine, bromine)
  - Aromatic or halogenated hydrocarbons
  - Corrosion inhibitors
- Do not use rinse aids or neutralisers.

#### **WARNING**

The use of incorrect temperatures may damage implants, instruments, caddies, transport racks and sterilization containers.

- Damaged implants and instruments may fail.
- Implants, instruments, caddies, transport racks and sterilization containers should only be exposed to a maximum temperature of 142 °C (288 °F).

### 7.6 Technical data

#### **ATTENTION**

Damage to implants and instruments through the use of incorrect operating resources.

- As regards operating resources, the requirements given in chapter 7.5 „Material resilience“ must also be observed.

### 7.6.1 Operating resources for pre-treatment of instruments

(according to chapter 7.1.2)

The suitability of the instruments for pre-treatment has been verified by an independent, accredited and officially recognized testing laboratory using the following cleaning agents:

Medium	Name	Manufacturer
Cleaning agent	Neodisher mediclean forte	Dr. Weigert GmbH & Co. KG, Hamburg

The procedure followed corresponds to the procedure described in chapter 7.1.2 „Pre-treatment in the Central Sterile Services Department“.

Properties of operating resources	
Cleaning agent <sup>1</sup>	<ul style="list-style-type: none"> <li>- Approved for cleaning instruments made of metal and plastic</li> <li>- Compatible with the materials of the instruments</li> <li>- Suitable for use with ultrasound</li> <li>- Free of alcohols and simple aldehydes, otherwise protein fixation</li> </ul>
Water	<ul style="list-style-type: none"> <li>- To max. 45 °C (113 °F)</li> <li>- Sterile or low-germ with max. 10 CFU/ml</li> <li>- Low in endotoxins at max. 0.25 EU/ml</li> <li>→ e.g. purified water or highly purified water</li> </ul>
Additional aids	<ul style="list-style-type: none"> <li>- Absorbent and low-linting single-use cloths</li> <li>- Soft brushes with varying diameters</li> <li>- Syringes with varying volumes (for details see chapter 7.1.8)</li> <li>- Cannulas of varying sizes (for details see chapter 7.1.8)</li> <li>- For the group allocation of instruments, observe the surgical instructions relevant to the implant system</li> <li>- Ultrasonic bath</li> </ul>

### 7.6.2 Operating resources for machine cleaning and disinfection of instruments

(according to chapter 7.1.4)

The suitability of the implants and instruments for machine cleaning and disinfection has been verified by an independent, accredited and officially recognized test lab using the following devices and cleaning agents:

Device	Method	Manufacturer
RDG G 7836 CD	Thermal disinfection	Miele & Cie. GmbH & Co., Gütersloh

Medium	Name	Manufacturer
Cleaning agent	Neodisher mediclean forte	Dr. Weigert GmbH & Co. KG, Hamburg

The procedure followed corresponds to the procedure described in chapter 7.1.4 „Machine cleaning and disinfection“.

Properties of operating resources	
Cleaning agent system	<ul style="list-style-type: none"> <li>- Approved for cleaning implants and instruments made of metal and plastic</li> <li>- Compatible with the materials of the implants and instruments</li> </ul>

<sup>1</sup> If you are using a cleaning and disinfection agent, for example for health and safety reasons, please bear in mind that it should be aldehyde-free (otherwise blood contaminants will become fixed), should offer proven efficacy (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE certification), is suitable for disinfection of the instruments and is compatible with them (see chapter 7.5 „Material resilience“). Please note that the disinfectants used in pre-treatment only serve to protect staff and cannot replace the subsequent – post-cleaning – disinfection step.

**Properties of operating resources**

Cleaning and disinfecting device	<ul style="list-style-type: none"> <li>- With proven efficacy (e.g. VAH/DGHM approval or FDA approval or CE certification as per EN ISO 15883)</li> <li>- Compliance with national regulations</li> <li>- Sufficient number of rinsing cycles</li> <li>- Use water that is sterile or low-germ and low in endotoxins <ul style="list-style-type: none"> <li>- Germs: Max. 10 CFU/ml</li> <li>- Endotoxins: max. 0.25 EU/ml → e.g. purified water or highly purified water</li> </ul> </li> <li>- Air for drying: filtered, oil-free, low-germ, low-particle</li> <li>- Suitable for thermal disinfection with A0 value of <math>\geq 3000</math></li> <li>- Suitable for implants and instruments made of metal and plastic</li> <li>- Always perform thermal disinfection</li> </ul>
----------------------------------	---

**7.6.3 Operating resources for the care of instruments**

(according to chapter 7.1.6)

The suitability of the procedure for the care of instruments has been verified by an independent, accredited and officially recognized test laboratory using the following instrument oil:

Medium	Name	Manufacturer
Instrument oil	Spezial Ölspray	Dr. Schumacher GmbH, Malsfeld

For the care of the instruments an oil for surgical instruments can be used to lubricate all moving parts of the instruments. This is at the discretion of the user.

**Properties of operating resources**

Instrument oil	- based on white paraffin oil
----------------	-------------------------------

** WARNING**

Damage and contamination of instruments due to use of unsuitable instrument oils.

- Residues can result in risks for the patient.
  - Damaged instruments may fail.
- Do not use oils with other additives such as silicones or corrosion inhibitors.  
→ Use instrument oil based on white paraffin oil.

**7.6.4 Operating resources for sterilization**

(according to chapter 7.2)

The suitability of the implants and instruments for sterilization has been verified by an independent, accredited and officially recognised test laboratory.

**Properties of operating resources**

Steam sterilizer	<ul style="list-style-type: none"> <li>- Design as per EN 13060 or EN 285 or ANSI/AAMI ST79 for USA: FDA cleared</li> <li>- Validation as per EN ISO 17665</li> <li>- Maximum sterilization temperature plus tolerance acc. to DIN EN ISO 17665: 138 °C (280 °F)</li> </ul>
------------------	---

<b>Properties of operating resources</b>	
Sterilization container	<ul style="list-style-type: none"><li>- Compliance with national regulations, including EN ISO / ANSI/AAMI/ISO 11607 / EN 868-8 for USA: FDA cleared</li><li>- Suitable for steam sterilization: To withstand temperatures up to at least 142 °C (288 °F), sufficiently steam permeable</li><li>- Protection of products and sterilization trays from mechanical damage</li></ul>
Wrapping material for containers without integrated filter system	<ul style="list-style-type: none"><li>- With proven efficacy (FDA approval or CE certification)</li><li>- With compliance with national regulations</li></ul>

## **8. NOTICE FOR THE SCOPE OF EU REGULATION (MDR) 2017/745, ANNEX I, CHAP. III z)**

All serious incidents that have occurred in connection with the product must be reported to the manufacturer and the competent authorities of the member state of the European Community in which the user and/or the patient is established, in accordance with EU Regulation (MDR) 2017/745.

## 9. SYMBOLS

The following table defines the symbols used for the product, labelling or the instruction for use of the system according to ISO 15223-1.

Symbol	Title	Description of symbol (acc. ISO 15223-1)
	Manufacturer	Indicates the medical device manufacturer.
	Date of manufacture	Indicates the date on which the medical device was manufactured.
	Medical device	Indicates that the item in question is a medical device.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
 <small>www.biedermann.com/eifu</small>	Consult electronic instructions for use	Indicates the need for the user to consult the instructions for use. <small>www.biedermann.com/eifu</small>
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.

The following symbol is used according to the FDA „Guidance for Industry – Alternative to Certain Prescription Device Labeling Requirements“.

Symbol	Title	Description
	Rx only	Caution: Federal law restricts this device to sale in the USA by or in the order of a physician.

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Biedermann Motech GmbH & Co. KG